



entermed

MICRONOMIC COMBI®

The MICRONOMIC 50® and MICRONOMIC PLUS® can be combined with connecting shelves and a long Marlan worktop to form the MICRONOMIC COMBI®

Approx dimensions of the MICRONOMIC COMBI®:

Height: 90 cm
Length: 160 cm
Depth: 60 cm
Weight: 120 kg (empty)

- Containers for waste and used instruments
- Handpiece holder with infrared sensor
- Wide choice of integrated functions
- Mobile on castors
- Available in 3 standard colours

- Compressed air system c/w 2 liquid 1 powder spray bottles
- Ear syringe with electronic temperature control
- Mirror preheater
- Hot air mirror heater
- Heated telescope holders
- Non-heated telescope holders
- Low voltage power source with switchable hook for Frenzel glasses and examination headlight
- Diagnostic microscope with 3-step magnification changer
- Full HD camera system with integral computer and USB connector in console
- 21" touch-screen HD monitor
- 100 Watt HF generator (mono and bipolar)
- X-ray viewer on a post



MICRONOMIC COMBI® with microscope, HD camera system and 100 W HF generator



Product code MICRONOMIC COMBI®: 980566