

Elektro-Mag M20-40

An intelligent, microprocessor controlled electrosurgical unit



Elektro-Mag M20-40 is an intelligent, microprocessor controlled electrosurgical unit. Based on years of clinic research, it is safe, reliable and powerful and can be used in the outpatients department as well as in the operating theatre. Elektro-Mag M20-40 is a versatile unit that has a wide spectrum of clinical applications (ENT, plastic surgery, proctology, dermatology, gynaecology, outpatient surgery) and guarantees fast tissue healing and quick patient recovery.

*A versatile unit that has a wide spectrum
of clinical applications*

Benefits at a glance

- 100 Watts maximum output power
- 9 presets
- 5 operating modes: spray coagulation; pure and blended cutting modes; bipolar and coagulation
- Finger and foot switch operation
- High level of patient safety



Technical specifications

	Maximal Output Power
Pure CUT	100W / 300 Ω
Blended CUT	100W / 300 Ω
Contact COAG	40W / 300 Ω
Spray COAG	30W / 1000 Ω
Bipolar COAG	20W / 75 Ω
Working frequency	400 kHz
Modulation frequency	30 kHz
Activation	Dual Footswitch and Pencil with switch
Dimensions (W x D x H)	300 mm x 240 mm x 115 mm
Weight	7.5 kg

Product code: M20-40T

Elektro-Mag M20-40 electro-surgical generator with: handle with 3m cable, double foot switch, 3m cable for neutral electrode and a set of monopolar electrodes and a bipolar forceps, patient plate, 3m bipolar cable, power cord and user manual