


# Endocontainer

Endoscope sterilisation container



The Endocontainer keeps endoscopes safe during the sterilisation process. The container accommodates two rigid endoscopes of up to 175 mm in length of various diameters. The soft silicone projections do not damage endoscope whilst holding them securely in position.

Robust endoscope sterilisation  
container

## Benefits at a glance

- Keeps endoscopes safe during the sterilisation process
- Robust stainless steel construction
- Sturdy lid with secure closing mechanism



Robust stainless steel construction

## Technical specifications

<b>Width</b>	80 mm
<b>Length</b>	290 mm
<b>Height</b>	50 mm

**Product code: cont. 29.8.5**

Endocontainer

Endoscopes sterilisation container