

Entermed Full HD camera

USB 3.0 Endoscopy camera



The Entermed Full HD Endoscopy Camera is a high definition camera designed for use with a PC and HD monitor. The camera comes with an easy-to-use software package, which allows images and videos to be saved to disk via the foot control provided. The software allows the brightness, contrast, hue, colour saturation and sharpness of the picture to be adjusted according to specific circumstances. Thumbnails make it possible to conveniently revert to saved images and videos during the session. The picture quality is remarkably crisp and clear.

The picture quality is
remarkably crisp and clear

Benefits at a glance

- Full HD camera USB 3.0
- CMOS progressive scan system
- Software includes anti-moiré filter
- User-friendly software
- Hands-free camera controls
- Automatic white balance
- Soakable camera head



Full HD camera USB 3.0

Technical specifications

Resolution	Full HD (1920 x 1080p / 30fps)
Cable length	2.5 m
Operating system	Windows 7 or later
Image device	1/2.8" CMOS SONY progressive scan system
Aspect ratio	16 : 9
Gain control	Automatic
Electronic shutter	1/50 ~ 1/100.000 sec
Supply voltage	USB 5.0V
Operating temperature	-20°C ~ 50°C
Storage temperature	-20°C ~ 50°C
Dimensions	diameter 35 mm, 120 mm length
Minimum illumination	1.0 Lux (F-12/20 IRE)

Product code: 840300

Entermed Full HD camera USB 3.0
with foot switch and imaging software

Product code: 830016

HD endoscopy adaptor with optical zoom 16 mm to 32 mm

Product code: 600493

24" monitor, full HD, touch screen