

ENTerview® III - 5500

LED cordless examination headlight



The ENTerview® III - 5500 is a rechargeable, cordless headlight which overcomes the need for connection to a fibre optic light source, separate battery pack or transformer, and hence increases freedom and mobility. The balanced design of the ENTerview® III - 5500 headlight permits comfortable use during longer periods of time. The light emitting diode (LED) provides ample light for a regular outpatient session. The ENTerview® III - 5500 is a versatile aid that is suitable for use in ENT, plastic surgery, gynaecology, maxillo facial surgery etc.

Automatically switches on when taken off its battery charger. Switches off when placed back on its battery charger after use

Benefits at a glance

- Cordless, rechargeable design
- Automatic on/off function
- Bright LED light
- Adjustable spot size
- Fully adjustable link between headband and lamp housing
- Adjustable headband with removable covering
- Multipurpose
- Light weight



Comfortable fit



Dimmer function



Adjustable spot size

Technical specifications

- Weight ca 305 grams
- Spot diameter (distance 25 cm) adjustable 25-100 mm
- Colour temperature 5500K
- Charging by specially designed battery charger
- Battery life sufficient for approx 4.5 hours use
- Luminous flux typical 235 lumen
- Collimating lens system
- Fresnel lens

Product code: 1030220-WAND-ENT

ENTERVIEW® III – 5500 cordless headlight
with wall-mounted battery charger

Product code: 1030220-ENT

ENTERVIEW® III – 5500 cordless headlight
with battery charger mounted on Entermed ENT treatment unit