

Flexible Naso-Pharyngo -Laryngoscope

Flexible endoscope



This superior flexible endoscope is manufactured in Germany to the highest of standards. It is highly suitable for both adult and paediatric use. The high resolution, crystal clear optics ensure sharp image quality. The flexible tip can be moved 150° using the scope's mechanism.

Highly suitable
for both adult and paediatric use

Benefits at a glance

- Crystal clear optics
- Precise motion for perfect control
- The flexible tip can be moved 150°
- Waterproof design
- Supplied in a carrying case with leak tester and Storz/Wolf adaptors



Leak tester

Technical specifications

ENTERMED naso-pharyngo-laryngo scope

- Working length 300mm
- Diameter of distal tip 2.8 mm or a 3.4 mm
- Bending 150° up, 150° down
- 70° angle of vision
- Depth of focus 5-50mm
- 18.000 pixels



Supplied in a carrying case

Product code:

700828E PAEDIATRIC

700800E ADULT

Flexible Naso - Pharyngo - Laryngoscope

Flexible endoscope